

pHDE-35S-Cas9-mCherry

Cat PP00435

Vector type

Plant Expression, CRISPR

Standard format

Vector backbone

pHDE

Modifications to backbone

mCherry expression cassette
35S-Cas9

Selectable markers

Hygromycin ; mCherry fluorescence

Bacterial Resistance(s)

Spectinomycin, 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DH5α

Copy number

Unknown

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY