

# N5m-GFP pZK538

Cat PP00476

---

## Vector type

Bacterial Expression, Plant Expression

## Standard format

## Total vector size (bp)

14169

## Vector backbone

pZK538

## Selectable markers

Phosphinothricin

## Bacterial Resistance(s)

Kanamycin, 50 µg/mL

## Growth Temperature

37°C

## Growth Strain(s)

DH5α

## Copy number

High Copy

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**