

# pEG302 22aa SunTag VP64 gRNA17 (FWA)

Cat *PP00487*

---

## Vector type

Plant Expression, CRISPR, Synthetic Biology

## Standard format

## Vector backbone

pEG302

## Selectable markers

Hygromycin

## Bacterial Resistance(s)

Kanamycin, 50 µg/mL

## Growth Temperature

37°C

## Growth Strain(s)

DH5α

## Copy number

Unknown

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**