

Erythrina cristagalli Lectin (ECL/ECA) - MagneZoom™

Cat PL00126

Size 2 mL/5 mL

Sugar Specificity

Gal β 4GlcNAc

Description

Affinity purified red algal lectins (ECL/ECA) are fixed on aminated paramagnetic beads with EDC and stabilized in buffers to maintain maximum activity. These beads are ideal for iP applications to pull down glycoproteins.

Abbreviation

ECL/ECA

Material Source

Coral tree seeds

Conjugate

MagneZoom™

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal β 4GlcNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY