

UBQ10:sXVE:ATG-3xHA-S11

Cat PP00756

Vector type

Plant Expression, Synthetic Biology

Standard format

Total vector size (bp)

14561

Vector backbone

UBQ10:sXVE:GFPc:Kan

Selectable markers

Kanamycin

Bacterial Resistance(s)

Kanamycin, 50 µg/mL

Growth Temperature

37°C

Growth Strain(s)

DH5α

Copy number

Low Copy

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY