

# Plant Total Protein Miniprep Kit (SDS-PAGE)

Cat PE00003

Size 30T

---

## Shelf Life

1 year

## Application

This product is used for rapid extraction of trace plant proteins for SDS-PAGE gel electrophoresis and Western blot analysis.

## Description

A wide variety of plant material can be processed, including fresh or freeze-dried material.

Simple to operate.

Effectively removes common contaminants from plant samples .

## Kit Composition

Solution A 30ml

Solution B 240ml

Solution C 1.5ml

Instruction manual 1 copy

## Storage

RT, where Solution A and Solution B are stored at -20°C.

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**