

Erythrina cristagalli Lectin (ECL/ECA) - Separopore® 6B SuperTrap™ Cartridge

Cat PL00132

Size 3 mL/6 mL/5 mL/10 mL

Sugar Specificity

Description

The crosslinked 6B Separopore® is properly coupled to the ECL/ECA, providing a fast-flowing cartridge that can capture the affinity-bound material, retain it during washing, and then elute the affinity-purified material.

Purity

High Purity Grade

Shelf Life

1 year

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY