

Euonymus europaeus Lectin (EEL/EEA) - FITC (Fluorescein)

Cat PL00140

Size 2 mg

Sugar Specificity

Gal α 3Gal

Description

The European EEL/EEA is labeled with fluorescein isothiocyanate (FITC) and combined with appropriate amount of fluorescent dye to provide the best staining properties for the lectin. FITC conjugates have been used in a variety of immunohistochemical and flow cytometry applications.

Abbreviation

EEL/EEA

Material Source

Spindle tree

Conjugate

FITC (Fluorescein)

Concentration

2 mg/mL

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal α 3Gal

Inhibiting or Eluting Sugar

Lactose

Divalent Ions

Ca $^{++}$, Zn $^{++}$

Mitogenic Activity

Yes

Format

Liquid

Storage

-20°C

Excitation

495 nm

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Euonymus europaeus Lectin (EEL/EEA) - FITC (Fluorescein)

Cat *PL00140*

Size *2 mg*

Emission

515 nm

Fluorescence

Green

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY
