

Euonymus europaeus Lectin (EEL/EEA) - Separopore® 4B

Cat PL00144

Size 2 mL/5 mL

Sugar Specificity

Gal α 3Gal

Description

Affinity purified European Euonymus lectins (EEL/EEA) are attached to Separopore®4B after CNBr coupling. The final product is stabilized in the buffer to maintain maximum activity. The gel can be filled with 1 ml or 5 ml.

Abbreviation

EEL/EEA

Material Source

Spindle tree seeds

Conjugate

Separopore® 4B

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

Gal α 3Gal

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY