

# polyphenol oxidase (PPO), ELISA Kit

Cat PEK00001

---

## Species

Plant

## Full Product Name

Plant polyphenol oxidase (PPO) ELISA Kit

## Product Gene Name

PPO elisa kit

## Assay Type

Quantitative Sandwich

## Detection Range

62.5U/ml-2000U/ml

## Sensitivity

10U/ml

## Intra-assay Precision

Intra-assay CV (%) is less than 15%.

## Inter-assay Precision

Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean × 100]. All CV% should be compared by concentration, not compared by OD values.

## Storage

Store all reagents at 2-8°C

## Protein Family

Polyphenol oxidase

## NCBI Accession #

EPS57764.1

## NCBI GI #

527183303

## NCBI Official Full Name

polyphenol oxidase, partial

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY