

Cytokinin CYT, ELISA Kit

Cat PEK00007

Species

Plant

Full Product Name

Plant Cytokinin CYT ELISA Kit

Product Gene Name

CYT elisa kit

Assay Type

Quantitative Sandwich

Detection Range

100 ng/mL-1.56 ng/mL

Sensitivity

Up to 0.5 ng/mL

Intra-assay Precision

Inter-assay Precision

Samples

Plant serum, plasma or cell culture supernatant

Storage

Store all reagents at 2-8°C.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY