

Wheat germ agglutinin, WGA, ELISA Kit

Cat PEK00013

Species

Plant

Full Product Name

Wheat germ agglutinin, WGA ELISA Kit

Product Gene Name

WGA elisa kit

Assay Type

Quantitative Sandwich

Detection Range

0.25ng/ml-8ng/ml

Sensitivity

0.1ng/ml

Intra-assay Precision

Intra-assay CV (%) is less than 15%.

Inter-assay Precision

Inter-assay CV (%) is less than 15%. $[CV(\%) = SD/mean \times 100]$. All CV% should be compared by concentration, not compared by OD values.

Storage

Store all reagents at 2-8°C

NCBI Accession

AAA34257.1

NCBI GI

170668

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY