

Sucrose Phosphate Synthase SPS, ELISA Kit

Cat PEK00015

Species

Plant

Full Product Name

Plant Sucrose Phosphate Synthase SPS ELISA Kit

Product Gene Name

SPS elisa kit

Assay Type

Quantitative Sandwich

Detection Range

20 ng/ml-0.312 ng/ml

Sensitivity

Up to 0.06 ng/ml

Intra-assay Precision

Inter-assay Precision

Specificity

No cross-reaction with other factors.

Samples

Plant serum, plasma or cell culture supernatant

Storage

Store all reagents at 2-8°C.

Protein Family

Probable sucrose-phosphate synthase

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY