

Vitamin A, VA, ELISA Kit

Cat PEK00028

Species

Plant

Full Product Name

Plant Vitamin A, VA ELISA Kit

Product Gene Name

VA elisa kit

Assay Type

Competitive

Detection Range

0.625 ng/ml -10 ng/ml.

Sensitivity

0.467 ng/ml.

Intra-assay Precision

Intra-assay Precision (Precision within an assay): CV%

Inter-assay Precision

Inter-assay Precision (Precision between assays): CV%

Specificity

This assay has high sensitivity and excellent specificity for detection of plant VA. No significant cross-reactivity or interference between plant VA and analogues was observed.

Samples

Plant tissue.

Storage

Unopened test kits should be stored at 2 to 8°C upon receipt. Please refer to pdf manual for further storage instructions.

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY