

# Wingless Type MMTV Integration Site Family, Member 3A (WNT3A), ELISA Kit

Cat PEK00040

---

## Species

Plant

## Full Product Name

Plant Wingless Type MMTV Integration Site Family, Member 3A (WNT3A) ELISA Kit

## Product Gene Name

WNT3A elisa kit

## Molecular Weight

42,905 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

EAW69829.1

## NCBI GI #

119590235

## NCBI GeneID

89780

## NCBI Official Full Name

wingless-type MMTV integration site family, member 3A

## NCBI Official Symbol

WNT3A

## NCBI Official Synonym Full Names

Wnt family member 3A

## NCBI Protein Information

protein Wnt-3a

## UniProt Gene Name

WNT3A

## UniProt Protein Name

Protein Wnt-3a

## UniProt Entry Name

WNT3A\_HUMAN

## UniProt Secondary Accession #

Q3SY79; Q3SY80; Q969P2

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Wingless Type MMTV Integration Site Family, Member 3A (WNT3A), ELISA Kit

Cat PEK00040

---

## UniProt Related Accession #

P56704

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**