

Alternative Oxidase (ATO), ELISA Kit

Cat PEK00097

Species

Plant

Full Product Name

Plant Alternative Oxidase (ATO) ELISA Kit

Product Gene Name

ATO elisa kit

Assay Type

Sandwich

Detection Range

20 ng/ml-0.312 ng/ml

Sensitivity

0.06 ng/ml.

Intra-assay Precision

Inter-assay Precision

Specificity

No cross-reaction with other factors.

Samples

Plant serum, plasma or Cell Culture Supernatant and organizations in the natural and recombinant ATO concentration

Storage

Store all reagents at 2-8°C.

Protein Family

Protein atonal

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY