

Prostate Specific Antigen (PSAg), ELISA Kit

Cat PEK00146

Species

Plant

Full Product Name

Plant Prostate Specific Antigen (PSAg) ELISA Kit

Product Gene Name

PSAg elisa kit

Molecular Weight

25,116 Da

Storage

Store all reagents at 2-8°C

Protein Family

Photosystem I reaction center

NCBI Accession

CAD54617.1

NCBI GI

24370944

NCBI GeneID

354

NCBI Official Full Name

prostate specific antigen

NCBI Official Symbol

KLK3

NCBI Official Synonym Full Names

kallikrein related peptidase 3

NCBI Official Synonym Symbols

APS; PSA; hK3; KLK2A1

NCBI Protein Information

prostate-specific antigen

UniProt Gene Name

KLK3

UniProt Synonym Gene Names

APS; PSA; Semin

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Prostate Specific Antigen (PSAg), ELISA Kit

Cat PEK00146

Prostate-specific antigen

UniProt Synonym Protein Names

Gamma-seminoprotein; Semin; Kallikrein-3; P-30 antigen; Semenogelase

UniProt Entry Name

KLK3_HUMAN

UniProt Primary Accession

Q8IXI4

UniProt Secondary Accession

Q15096; Q16272; Q86TG8; Q8IXI4; C9JXH3; G3V0H4; G3XAE3

UniProt Related Accession

P07288

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY