

# Heparan Sulfate Proteoglycan (HSPG), ELISA Kit

Cat PEK00325

---

## Species

Plant

## Full Product Name

Plant Heparan Sulfate Proteoglycan (HSPG) ELISA Kit

## Product Gene Name

HSPG elisa kit

## Molecular Weight

468,830 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

AAA52700.1

## NCBI GI #

184427

## NCBI GeneID

3339

## NCBI Official Full Name

heparan sulfate proteoglycan

## NCBI Official Symbol

HSPG2

## NCBI Official Synonym Full Names

heparan sulfate proteoglycan 2

## NCBI Official Synonym Symbols

PLC; SJA; SJS; HSPG; SJS1; PRCAN

## NCBI Protein Information

basement membrane-specific heparan sulfate proteoglycan core protein

## UniProt Gene Name

HSPG2

## UniProt Synonym Gene Names

HSPG; PLC

## UniProt Protein Name

Basement membrane-specific heparan sulfate proteoglycan core protein

## UniProt Synonym Protein Names

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Heparan Sulfate Proteoglycan (HSPG), ELISA Kit

Cat PEK00325



---

Perlecan; PLC

## UniProt Entry Name

PGBM\_HUMAN

## UniProt Secondary Accession #

Q16287; Q5SZI3; Q9H3V5

## UniProt Related Accession #

P98160

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**