

Adenosine A2B Receptor (ADORA2B), ELISA Kit

Cat PEK00332

Species

Plant

Full Product Name

Plant Adenosine A2B Receptor (ADORA2B) ELISA Kit

Product Gene Name

ADORA2B elisa kit

Molecular Weight

36,333 Da

Storage

Store all reagents at 2-8°C

Protein Family

Adenosine receptor

NCBI Accession

AAA51598.1

NCBI GI

178150

NCBI GeneID

136

NCBI Official Full Name

adenosine A2b receptor

NCBI Official Symbol

ADORA2B

NCBI Official Synonym Full Names

adenosine A2b receptor

NCBI Official Synonym Symbols

ADORA2

NCBI Protein Information

adenosine receptor A2b

UniProt Gene Name

ADORA2B

UniProt Protein Name

Adenosine receptor A2b

UniProt Entry Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Adenosine A2B Receptor (ADORA2B), ELISA Kit

Cat PEK00332

AA2BR_HUMAN

UniProt Related Accession

P29275

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY