

Synaptonemal Complex Protein 3 (SYCP3), ELISA Kit

Cat PEK00340

Species

Plant

Full Product Name

Plant Synaptonemal Complex Protein 3 (SYCP3) ELISA Kit

Product Gene Name

SYCP3 elisa kit

Molecular Weight

27,729 Da

Storage

Store all reagents at 2-8°C

Protein Family

Synaptonemal complex protein

NCBI Accession

NP_001171420.1

NCBI GI

295789052

NCBI GeneID

50511

NCBI Official Full Name

synaptonemal complex protein 3

NCBI Official Symbol

SYCP3

NCBI Official Synonym Full Names

synaptonemal complex protein 3

NCBI Official Synonym Symbols

COR1; SCP3; SPGF4; RPRGL4

NCBI Protein Information

synaptonemal complex protein 3

UniProt Gene Name

SYCP3

UniProt Synonym Gene Names

SCP3; SCP-3

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Synaptonemal Complex Protein 3 (SYCP3), ELISA Kit

Cat PEK00340



Synaptonemal complex protein 3

UniProt Entry Name

SYCP3_HUMAN

UniProt Related Accession

Q8IZU3

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY