

# Proteinase 3 (PR3), ELISA Kit

Cat PEK00345

## Species

Plant

## Full Product Name

Plant Proteinase 3 (PR3) ELISA Kit

## Product Gene Name

PR3 elisa kit

## Molecular Weight

27,807 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

CAA39203.1

## NCBI GI #

1335280

## NCBI GeneID

5657

## NCBI Official Full Name

proteinase 3

## NCBI Official Symbol

PRTN3

## NCBI Official Synonym Full Names

proteinase 3

## NCBI Official Synonym Symbols

MBN; MBT; NP4; P29; PR3; ACPA; AGP7; NP-4; PR-3; CANCA; C-ANCA

## NCBI Protein Information

myeloblastin

## UniProt Gene Name

PRTN3

## UniProt Synonym Gene Names

MBN; PR-3; PR3; NP-4

## UniProt Protein Name

Myeloblastin

## UniProt Synonym Protein Names

AGP7; C-ANCA antigen; Leukocyte proteinase 3; PR-3; PR3; Neutrophil proteinase 4; NP-4; P29; Wegener autoantigen

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Proteinase 3 (PR3), ELISA Kit

Cat PEK00345

---

## UniProt Entry Name

PRTN3\_HUMAN

## UniProt Primary Accession #

P24158

## UniProt Secondary Accession #

P15637; P18078; Q4VB08; Q4VB09; Q6LBM7; Q6LBN2; Q9UD25; Q9UQD8

## UniProt Related Accession #

P24158

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**