

Triglyceride Lipase (TGL), ELISA Kit

Cat PEK00682

Species

Plant

Full Product Name

Plant Triglyceride Lipase (TGL) ELISA Kit

Product Gene Name

TGL elisa kit

Molecular Weight

51,157 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

AAA36740.1

NCBI GI

339597

NCBI GenelD

5406

NCBI Official Full Name

triglyceride lipase

NCBI Official Symbol

PNLIP

NCBI Official Synonym Full Names

pancreatic lipase

NCBI Official Synonym Symbols

PL; PTL; PNLIPD

NCBI Protein Information

pancreatic triacylglycerol lipase

UniProt Gene Name

PNLIP

UniProt Synonym Gene Names

PL; PTL; Pancreatic lipase

UniProt Protein Name

Pancreatic triacylglycerol lipase

UniProt Entry Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Triglyceride Lipase (TGL), ELISA Kit

Cat PEK00682

LIPP_HUMAN

UniProt Secondary Accession

Q5VSQ2

UniProt Related Accession

P16233

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY