

# High Mobility Group Box 1 Protein (HMGB1), ELISA Kit

Cat PEK00690

---

## Species

Plant

## Full Product Name

Plant High Mobility Group Box 1 Protein (HMGB1) ELISA Kit

## Product Gene Name

HMGB1 elisa kit

## Molecular Weight

24,894 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

High mobility group protein

## NCBI Accession #

ABM47301.1

## NCBI GI #

121073803

## NCBI GeneID

3146

## NCBI Official Full Name

high mobility group box 1 protein

## NCBI Official Symbol

HMGB1

## NCBI Official Synonym Full Names

high mobility group box 1

## NCBI Official Synonym Symbols

HMG1; HMG3; HMG-1; SBP-1

## NCBI Protein Information

high mobility group protein B1

## UniProt Gene Name

HMGB1

## UniProt Synonym Gene Names

HMG1; HMG-1

## UniProt Protein Name

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

---

# High Mobility Group Box 1 Protein (HMGB1), ELISA Kit

Cat PEK00690

---

High mobility group protein B1

## UniProt Synonym Protein Names

High mobility group protein 1; HMG-1

## UniProt Entry Name

HMGB1\_HUMAN

## UniProt Secondary Accession #

Q14321; Q5T7C3; Q6IBE1; A5D8W9

## UniProt Related Accession #

P09429

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**