

Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1), ELISA Kit

Cat PEK01004

Species

Plant

Full Product Name

Plant Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1) ELISA Kit

Product Gene Name

ADCYAP1 elisa kit

Molecular Weight

18,835 Da

Storage

Store all reagents at 2-8°C

Protein Family

Glucagon family neuropeptides

NCBI Accession

EAX01706.1

NCBI GI

119622111

NCBI GeneID

116

NCBI Official Full Name

adenylate cyclase activating polypeptide 1 (pituitary)

NCBI Official Symbol

ADCYAP1

NCBI Official Synonym Full Names

adenylate cyclase activating polypeptide 1

NCBI Official Synonym Symbols

PACAP

NCBI Protein Information

pituitary adenylate cyclase-activating polypeptide

UniProt Gene Name

ADCYAP1

UniProt Synonym Gene Names

PACAP; PACAP-27; PACAP27; PACAP-38; PACAP38

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Adenylate Cyclase Activating Polypeptide 1, Pituitary (ADCYAP1), ELISA Kit

Cat PEK01004

UniProt Protein Name

Pituitary adenylate cyclase-activating polypeptide

UniProt Synonym Protein Names

PACAP-related peptide
Alternative name(s): PRP-48

UniProt Entry Name

PACA_HUMAN

UniProt Secondary Accession

Q52LQ0; B2R7N4

UniProt Related Accession

P18509

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY