

# WNT Inhibitory Factor 1 (WIF1), ELISA Kit

Cat PEK01020

---

## Species

Plant

## Full Product Name

Plant WNT Inhibitory Factor 1 (WIF1) ELISA Kit

## Product Gene Name

WIF1 elisa kit

## Molecular Weight

41,590 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Wnt inhibitory factor

## NCBI Accession #

NP\_036045.1

## NCBI GI #

6755999

## NCBI GeneID

24117

## NCBI Official Full Name

wnt inhibitory factor 1

## NCBI Official Symbol

Wif1

## NCBI Official Synonym Full Names

Wnt inhibitory factor 1

## NCBI Official Synonym Symbols

WIF-1; AW107799

## NCBI Protein Information

wnt inhibitory factor 1

## UniProt Gene Name

Wif1

## UniProt Synonym Gene Names

WIF-1

## UniProt Protein Name

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

# WNT Inhibitory Factor 1 (WIF1), ELISA Kit

Cat PEK01020

---

Wnt inhibitory factor 1

## UniProt Entry Name

WIF1\_MOUSE

## UniProt Related Accession #

Q9WUA1

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**