

Activated Leukocyte Cell Adhesion Molecule (ALCAM), ELISA Kit

Cat PEK01041

Species

Plant

Full Product Name

Plant Activated Leukocyte Cell Adhesion Molecule (ALCAM) ELISA Kit

Product Gene Name

ALCAM elisa kit

Molecular Weight

33,878 Da

Storage

Store all reagents at 2-8°C

Protein Family

CD166 antigen

NCBI Accession

AAI37097.1

NCBI GI

187951595

NCBI GeneID

214

NCBI Official Full Name

Activated leukocyte cell adhesion molecule

NCBI Official Symbol

ALCAM

NCBI Official Synonym Full Names

activated leukocyte cell adhesion molecule

NCBI Official Synonym Symbols

MEMD; CD166

NCBI Protein Information

CD166 antigen

UniProt Gene Name

ALCAM

UniProt Protein Name

CD166 antigen

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Activated Leukocyte Cell Adhesion Molecule (ALCAM), ELISA Kit

Cat PEK01041

Activated leukocyte cell adhesion molecule; CD_antigen: CD166

UniProt Entry Name

CD166_HUMAN

UniProt Secondary Accession

O60892; Q1HGM8; Q1HGM9; Q6PEY4; Q6ZS95; B2RNS3; B4DTU0

UniProt Related Accession

Q13740

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY