

DNA Ligase 4 (LIG4), ELISA Kit

Cat PEK01054

Species

Plant

Full Product Name

Plant DNA Ligase 4 (LIG4) ELISA Kit

Product Gene Name

LIG4 elisa kit

Molecular Weight

103,971 Da

Storage

Store all reagents at 2-8°C

Protein Family

DNA ligase

NCBI Accession

NP_996820.1

NCBI GI

46255052

NCBI GeneID

3981

NCBI Official Full Name

DNA ligase 4

NCBI Official Symbol

LIG4

NCBI Official Synonym Full Names

DNA ligase 4

NCBI Official Synonym Symbols

LIG4S

NCBI Protein Information

DNA ligase 4

UniProt Gene Name

LIG4

UniProt Protein Name

DNA ligase 4

UniProt Synonym Protein Names

DNA ligase IV; Polydeoxyribonucleotide synthase [ATP] 4

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

DNA Ligase 4 (LIG4), ELISA Kit

Cat PEK01054

UniProt Entry Name

DNLI4_HUMAN

UniProt Secondary Accession

Q8IY66; Q8TEU5

UniProt Related Accession

P49917

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY