

Phospholipid Transfer Protein (PLTP), ELISA Kit

Cat PEK01065

Species

Plant

Full Product Name

Plant Phospholipid Transfer Protein (PLTP) ELISA Kit

Product Gene Name

PLTP elisa kit

Molecular Weight

44,847 Da

Storage

Store all reagents at 2-8°C

Protein Family

Phospholipid transfer protein

NCBI Accession

AAA36443.1

NCBI GI

468326

NCBI GeneID

5360

NCBI Official Full Name

phospholipid transfer protein

NCBI Official Symbol

PLTP

NCBI Official Synonym Full Names

phospholipid transfer protein

NCBI Official Synonym Symbols

BPIFE; HDLCQ9

NCBI Protein Information

phospholipid transfer protein

UniProt Gene Name

PLTP

UniProt Protein Name

Phospholipid transfer protein

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Phospholipid Transfer Protein (PLTP), ELISA Kit

Cat PEK01065

Lipid transfer protein II

UniProt Entry Name

PLTP_HUMAN

UniProt Secondary Accession

Q8WTT1; Q9BR07; Q9BSH8; A8K006; B4DDD5; B4DRB4; E1P5N8; E7EV16

UniProt Related Accession

P55058

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY