

Dopamine Receptor D2 (DRD2), ELISA Kit

Cat PEK01075

Species

Plant

Full Product Name

Plant Dopamine Receptor D2 (DRD2) ELISA Kit

Product Gene Name

DRD2 elisa kit

Molecular Weight

50,847 Da

Storage

Store all reagents at 2-8°C

Protein Family

D(2) dopamine receptor

NCBI Accession

CAB56463.1

NCBI GI

5921480

NCBI GeneID

1813

NCBI Official Full Name

dopamine receptor D2

NCBI Official Symbol

DRD2

NCBI Official Synonym Full Names

dopamine receptor D2

NCBI Official Synonym Symbols

D2R; D2DR

NCBI Protein Information

D(2) dopamine receptor

UniProt Gene Name

DRD2

UniProt Protein Name

D(2) dopamine receptor

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Dopamine Receptor D2 (DRD2), ELISA Kit

Cat PEK01075

Dopamine D2 receptor

UniProt Entry Name

DRD2_HUMAN

UniProt Primary Accession

P14416

UniProt Secondary Accession

Q9NZR3; Q9UPA9

UniProt Related Accession

P14416

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY