

# Programmed Cell Death Protein 1 (PDCD1), ELISA Kit

Cat PEK01086

## Species

Plant

## Full Product Name

Plant Programmed Cell Death Protein 1 (PDCD1) ELISA Kit

## Product Gene Name

PDCD1 elisa kit

## Molecular Weight

31,647 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Programmed cell death protein

## NCBI Accession #

NP\_005009.2

## NCBI GI #

167857792

## NCBI GeneID

5133

## NCBI Official Full Name

programmed cell death protein 1

## NCBI Official Symbol

PDCD1

## NCBI Official Synonym Full Names

programmed cell death 1

## NCBI Official Synonym Symbols

PD1; PD-1; CD279; SLEB2; hPD-1; hPD-I; hSLE1

## NCBI Protein Information

programmed cell death protein 1

## UniProt Gene Name

PDCD1

## UniProt Synonym Gene Names

PD1; Protein PD-1; hPD-1

## UniProt Protein Name

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Programmed Cell Death Protein 1 (PDCD1), ELISA Kit

Cat PEK01086

---

Programmed cell death protein 1

## UniProt Synonym Protein Names

CD\_antigen: CD279

## UniProt Entry Name

PDCD1\_HUMAN

## UniProt Secondary Accession #

O00517; Q8IX89

## UniProt Related Accession #

Q15116

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**