

Fibronectin Leucine Rich Transmembrane Protein 2 (FLRT2), ELISA Kit

Cat PEK01093

Species

Plant

Full Product Name

Plant Fibronectin Leucine Rich Transmembrane Protein 2 (FLRT2) ELISA Kit

Product Gene Name

FLRT2 elisa kit

Molecular Weight

74,049 Da

Storage

Store all reagents at 2-8°C

Protein Family

Leucine-rich repeat transmembrane protein

NCBI Accession

AAI30291.1

NCBI GI

120660314

NCBI GeneID

23768

NCBI Official Full Name

Fibronectin leucine rich transmembrane protein 2

NCBI Official Symbol

FLRT2

NCBI Official Synonym Full Names

fibronectin leucine rich transmembrane protein 2

NCBI Protein Information

leucine-rich repeat transmembrane protein FLRT2

UniProt Gene Name

FLRT2

UniProt Synonym Gene Names

KIAA0405

UniProt Protein Name

Leucine-rich repeat transmembrane protein FLRT2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Fibronectin Leucine Rich Transmembrane Protein 2 (FLRT2), ELISA Kit

Cat PEK01093



UniProt Synonym Protein Names

Fibronectin-like domain-containing leucine-rich transmembrane protein 2

UniProt Entry Name

FLRT2_HUMAN

UniProt Secondary Accession

A0AV84; B7ZLP3

UniProt Related Accession

O43155

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY