

# Plasminogen (PLG), ELISA Kit

Cat PEK01117

## Species

Plant

## Full Product Name

Plant Plasminogen (PLG) ELISA Kit

## Product Gene Name

PLG elisa kit

## Molecular Weight

90,569 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Plasminogen

## NCBI Accession #

AAA60113.1

## NCBI GI #

387026

## NCBI Official Full Name

plasminogen

## UniProt Gene Name

PLG

## UniProt Protein Name

Plasminogen

## UniProt Entry Name

PLMN\_HUMAN

## UniProt Secondary Accession #

Q15146; Q5TEH4; Q6PA00

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**