

Phospholipid Scramblase 2 (PLSCR2), ELISA Kit

Cat PEK01118

Species

Plant

Full Product Name

Plant Phospholipid Scramblase 2 (PLSCR2) ELISA Kit

Product Gene Name

PLSCR2 elisa kit

Molecular Weight

34,052 Da

Storage

Store all reagents at 2-8°C

Protein Family

Phospholipid scramblase

NCBI Accession

NP_032906.2

NCBI GI

31982115

NCBI GeneID

18828

NCBI Official Full Name

phospholipid scramblase 2

NCBI Official Symbol

Plscr2

NCBI Official Synonym Full Names

phospholipid scramblase 2

NCBI Protein Information

phospholipid scramblase 2

UniProt Gene Name

Plscr2

UniProt Synonym Gene Names

PL scramblase 2

UniProt Protein Name

Phospholipid scramblase 2

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Phospholipid Scramblase 2 (PLSCR2), ELISA Kit

Cat PEK01118

Ca(2+)-dependent phospholipid scramblase 2

UniProt Entry Name

PLS2_MOUSE

UniProt Secondary Accession

O70233; Q921P6

UniProt Related Accession

Q9DCW2

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY