

# X-Linked Inhibitor of Apoptosis Protein (XIAP), ELISA Kit

Cat PEK01379

---

## Species

Plant

## Full Product Name

Plant X-Linked Inhibitor of Apoptosis Protein (XIAP) ELISA Kit

## Product Gene Name

XIAP elisa kit

## Molecular Weight

56,685 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

E3 ubiquitin-protein ligase

## NCBI Accession #

AAC50373.1

## NCBI GI #

1184320

## NCBI Official Full Name

X-linked inhibitor of apoptosis protein

## UniProt Gene Name

XIAP

## UniProt Synonym Gene Names

API3; BIRC4; IAP3; ILP; hILP; IAP-3; hIAP-3; hIAP3; X-linked IAP

## UniProt Protein Name

E3 ubiquitin-protein ligase XIAP

## UniProt Synonym Protein Names

Baculoviral IAP repeat-containing protein 4; IAP-like protein; ILP; hILP; Inhibitor of apoptosis protein 3; IAP-3; hIAP-3; hIAP3; X-linked inhibitor of apoptosis protein; X-linked IAP

## UniProt Entry Name

XIAP\_HUMAN

## UniProt Secondary Accession #

Q9NQ14; D3DTF2

---

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY