

Glycine max Lectin (SBA) - Separopore® 4B

Cat PL00162

Size 2 mL/5 mL

Sugar Specificity

α > β GalNAc

Description

Affinity purified Glycine max lectin (SBA) is attached to Separopore®4B after CNBr coupling. The final product is stabilized in the buffer to maintain maximum activity. The gel can be filled with 1 ml or 5 ml.

Abbreviation

SBA

Material Source

Soybean seeds

Conjugate

Separopore® 4B

Purity

High Purity Grade

Shelf Life

1 year

Blood Group Specificity

Non-specific

Sugar Specificity

α > β GalNAc

Mitogenic Activity

No

Format

Liquid

Storage

2-8°C

Hazardous Shipping

Non-hazardous

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY