

Peptidyl Arginine Deiminase Type IV (PADI4), ELISA Kit

Cat PEK01708

Species

Plant

Full Product Name

Plant Peptidyl Arginine Deiminase Type IV (PADI4) ELISA Kit

Product Gene Name

PADI4 elisa kit

Molecular Weight

10,961 Da

Storage

Store all reagents at 2-8°C

Protein Family

Protein-arginine deiminase

NCBI Accession

CAG30581.2

NCBI GI

50299840

NCBI GeneID

23569

NCBI Official Full Name

peptidyl arginine deiminase type IV, partial

NCBI Official Symbol

PADI4

NCBI Official Synonym Full Names

peptidyl arginine deiminase 4

NCBI Official Synonym Symbols

PAD; PAD4; PDI4; PDI5; PADI5

NCBI Protein Information

protein-arginine deiminase type-4

UniProt Gene Name

PADI4

UniProt Protein Name

Peptidyl arginine deiminase type IV

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Peptidyl Arginine Deiminase Type IV (PADI4), ELISA Kit

Cat PEK01708

Peptidyl arginine deiminase type IV

UniProt Entry Name

Q6EVJ7_HUMAN

UniProt Primary Accession

Q6EVJ7

UniProt Related Accession

Q9UM07

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY