

# Interleukin 17 Receptor A (IL17RA), ELISA Kit

Cat PEK01710

---

## Species

Plant

## Full Product Name

Plant Interleukin 17 Receptor A (IL17RA) ELISA Kit

## Product Gene Name

IL17RA elisa kit

## Molecular Weight

92,153 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Interleukin-17 receptor

## NCBI Accession #

AAH11624.1

## NCBI GI #

15079611

## NCBI GeneID

23765

## NCBI Official Full Name

Interleukin 17 receptor A

## NCBI Official Symbol

IL17RA

## NCBI Official Synonym Full Names

interleukin 17 receptor A

## NCBI Official Synonym Symbols

CD217; IL17R; CANDF5; CDw217; IL-17RA; hIL-17R

## NCBI Protein Information

interleukin-17 receptor A

## UniProt Gene Name

IL17RA

## UniProt Synonym Gene Names

IL17R; IL-17 receptor A; IL-17RA

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Interleukin 17 Receptor A (IL17RA), ELISA Kit

Cat PEK01710

---

Interleukin-17 receptor A

## UniProt Synonym Protein Names

CDw217; CD\_antigen: CD217

## UniProt Entry Name

I17RA\_HUMAN

## UniProt Secondary Accession #

O43844; Q20WK1

## UniProt Related Accession #

Q96F46

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**