

# Dual Specificity Phosphatase 6 (DUSP6), ELISA Kit

Cat PEK01721

---

## Species

Plant

## Full Product Name

Plant Dual Specificity Phosphatase 6 (DUSP6) ELISA Kit

## Product Gene Name

DUSP6 elisa kit

## Molecular Weight

26,468 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Dual specificity protein phosphatase

## NCBI Accession #

AAP35541.1

## NCBI GI #

30582629

## NCBI GeneID

1848

## NCBI Official Full Name

dual specificity phosphatase 6

## NCBI Official Symbol

DUSP6

## NCBI Official Synonym Full Names

dual specificity phosphatase 6

## NCBI Official Synonym Symbols

HH19; MKP3; PYST1

## NCBI Protein Information

dual specificity protein phosphatase 6

## UniProt Gene Name

DUSP6

## UniProt Synonym Gene Names

MKP3; PYST1; MAP kinase phosphatase 3; MKP-3

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Dual Specificity Phosphatase 6 (DUSP6), ELISA Kit

Cat PEK01721

---

Dual specificity protein phosphatase 6

## UniProt Synonym Protein Names

Dual specificity protein phosphatase PYST1; Mitogen-activated protein kinase phosphatase 3; MAP kinase phosphatase 3; MKP-3

## UniProt Entry Name

DUS6\_HUMAN

## UniProt Secondary Accession #

O75109; Q53Y75; Q9BSH6

## UniProt Related Accession #

Q16828

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**