

# Aquaporin 4 (AQP-4), ELISA Kit

Cat PEK01726

## Species

Plant

## Full Product Name

Plant Aquaporin 4 (AQP-4) ELISA Kit

## Product Gene Name

AQP-4 elisa kit

## Molecular Weight

32,300 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

AAH22286.1

## NCBI GI #

18490380

## NCBI GeneID

361

## NCBI Official Full Name

Aquaporin 4

## NCBI Official Symbol

AQP4

## NCBI Official Synonym Full Names

aquaporin 4

## NCBI Official Synonym Symbols

MIWC; WCH4

## NCBI Protein Information

aquaporin-4

## UniProt Gene Name

AQP4

## UniProt Synonym Gene Names

AQP-4; MIWC

## UniProt Protein Name

Aquaporin-4

## UniProt Synonym Protein Names

Mercurial-insensitive water channel; MIWC; WCH4

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# Aquaporin 4 (AQP-4), ELISA Kit

Cat PEK01726

---

## UniProt Entry Name

AQP4\_HUMAN

## UniProt Secondary Accession #

P78564

## UniProt Related Accession #

P55087

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**