

Interferon Alpha 2A (IFNA2A), ELISA Kit

Cat PEK01745

Species

Plant

Full Product Name

Plant Interferon Alpha 2A (IFNA2A) ELISA Kit

Product Gene Name

IFNA2A elisa kit

Molecular Weight

21,550 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

NP_000596.2

NCBI GI

11067751

NCBI GeneID

3440

NCBI Official Full Name

interferon alpha-2

NCBI Official Symbol

IFNA2

NCBI Official Synonym Full Names

interferon alpha 2

NCBI Official Synonym Symbols

IFNA; INFA2; IFNA2B; IFN-alphaA

NCBI Protein Information

interferon alpha-2

UniProt Gene Name

IFNA2

UniProt Synonym Gene Names

IFNA2A; IFNA2B; IFNA2C; IFN-alpha-2; LeIF A

UniProt Protein Name

Interferon alpha-2

UniProt Synonym Protein Names

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Interferon Alpha 2A (IFNA2A), ELISA Kit

Cat PEK01745

Interferon alpha-A; LeIF A

UniProt Entry Name

IFNA2_HUMAN

UniProt Secondary Accession

P01564; Q14606; Q6DJX8; Q96KI6; H2DF54; H2DF55

UniProt Related Accession

P01563

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY