

# GM2 Ganglioside Activator Protein (GM2AP), ELISA Kit

Cat PEK01749

---

## Species

Plant

## Full Product Name

Plant GM2 Ganglioside Activator Protein (GM2AP) ELISA Kit

## Product Gene Name

GM2AP elisa kit

## Molecular Weight

20,610 Da

## Storage

Store all reagents at 2-8°C

## NCBI Accession #

NP\_001270491.1

## NCBI GI #

545686971

## NCBI GeneID

102143555

## NCBI Official Full Name

ganglioside GM2 activator

## NCBI Official Symbol

GM2A

## NCBI Official Synonym Symbols

SAP-3; GM2-AP

## NCBI Protein Information

ganglioside GM2 activator

## UniProt Gene Name

GM2A

## UniProt Synonym Gene Names

SAP-3

## UniProt Protein Name

Ganglioside GM2 activator

## UniProt Synonym Protein Names

Cerebroside sulfate activator protein; GM2-AP; Sphingolipid activator protein 3; SAP-3

## UniProt Entry Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# GM2 Ganglioside Activator Protein (GM2AP), ELISA Kit



Cat PEK01749

---

SAP3\_MACFA

## UniProt Secondary Accession #

Q60HF5

## UniProt Related Accession #

Q8HXX6

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**