

Regulated On Activation In Normal T-Cell Expressed and Secreted (RANTES), ELISA Kit

Cat PEK01762

Species

Plant

Full Product Name

Plant Regulated On Activation In Normal T-Cell Expressed and Secreted (RANTES) ELISA Kit

Product Gene Name

RANTES elisa kit

Molecular Weight

9,990 Da

Storage

Store all reagents at 2-8°C

NCBI Accession

AAC03541.1

NCBI GI

2905632

NCBI GeneID

6352

NCBI Official Full Name

RANTES

NCBI Official Symbol

CCL5

NCBI Official Synonym Full Names

C-C motif chemokine ligand 5

NCBI Official Synonym Symbols

SISd; eoCP; SCYA5; RANTES; TCP228; D17S136E; SIS-delta

NCBI Protein Information

C-C motif chemokine 5

UniProt Gene Name

CCL5

UniProt Synonym Gene Names

D17S136E; SCYA5; TCP228

UniProt Protein Name

C-C motif chemokine 5

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Regulated On Activation In Normal T-Cell Expressed and Secreted (RANTES), ELISA Kit

Cat PEK01762

UniProt Synonym Protein Names

EoCP; Eosinophil chemotactic cytokine; SIS-delta; Small-inducible cytokine A5; T cell-specific protein P228; TCP228

UniProt Entry Name

CCL5_HUMAN

UniProt Secondary Accession

O43646; Q0QVW8; Q4ZGJ1; Q9NYA2; Q9UBG2; Q9UC99

UniProt Related Accession

P13501

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY