

Paired Box Protein Pax-6 (PAX6), ELISA Kit

Cat PEK01765

Species

Plant

Full Product Name

Plant Paired Box Protein Pax-6 (PAX6) ELISA Kit

Product Gene Name

PAX6 elisa kit

Molecular Weight

49,961 Da

Storage

Store all reagents at 2-8°C

Protein Family

Paired box protein

NCBI Accession

NP_571379.1

NCBI GI

18859209

NCBI GeneID

30567

NCBI Official Full Name

paired box protein Pax-6

NCBI Official Symbol

pax6a

NCBI Official Synonym Full Names

paired box 6a

NCBI Official Synonym Symbols

pax6; cb280; pax-a; Pax6.1; paxzfa; PAXZF-A; fc20e07; zfpax-6a; zfpax-6b; pax[zf-a]; wu:fc20e07

NCBI Protein Information

paired box protein Pax-6

UniProt Gene Name

pax6a

UniProt Synonym Gene Names

pax[zf-a]; paxzf-a

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Paired Box Protein Pax-6 (PAX6), ELISA Kit

Cat PEK01765

Paired box protein Pax-6

UniProt Synonym Protein Names

Pax[Zf-a]

UniProt Entry Name

PAX6_DANRE

UniProt Secondary Accession

A2AUN5

UniProt Related Accession

P26630

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY