

# Eosinophil-Associated Ribonuclease A Family Member 3 (EAR3), ELISA Kit

Cat PEK02369

---

## Species

Plant

## Full Product Name

Plant Eosinophil-Associated Ribonuclease A Family Member 3 (EAR3) ELISA Kit

## Product Gene Name

EAR3 elisa kit

## Molecular Weight

17,745 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Eosinophil cationic-type ribonuclease

## NCBI Accession #

AAI53205.1

## NCBI GI #

157169982

## NCBI GenelD

53876

## NCBI Official Full Name

Eosinophil-associated, ribonuclease A family, member 3

## NCBI Official Symbol

Ear3

## NCBI Official Synonym Full Names

eosinophil-associated, ribonuclease A family, member 3

## NCBI Official Synonym Symbols

mR3; Rnase3

## NCBI Protein Information

eosinophil cationic-type ribonuclease 3

## UniProt Gene Name

Ear3

## UniProt Synonym Gene Names

Rnase3

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# Eosinophil-Associated Ribonuclease A Family Member 3 (EAR3), ELISA Kit

Cat PEK02369

---

## UniProt Protein Name

Eosinophil cationic-type ribonuclease 3

## UniProt Synonym Protein Names

MR-3

## UniProt Entry Name

ECP3\_MOUSE

## UniProt Related Accession #

O35290

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**