

Delta Like Protein 4 (DLL4), ELISA Kit

Cat PEK02373

Species

Plant

Full Product Name

Plant Delta Like Protein 4 (DLL4) ELISA Kit

Product Gene Name

DLL4 elisa kit

Molecular Weight

74,605 Da

Storage

Store all reagents at 2-8°C

Protein Family

Delta-like protein

NCBI Accession

NP_061947.1

NCBI GI

9506545

NCBI GeneID

54567

NCBI Official Full Name

delta-like protein 4

NCBI Official Symbol

DLL4

NCBI Official Synonym Full Names

delta like canonical Notch ligand 4

NCBI Official Synonym Symbols

AOS6; hdelta2

NCBI Protein Information

delta-like protein 4

UniProt Gene Name

DLL4

UniProt Synonym Gene Names

Delta4

UniProt Protein Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Delta Like Protein 4 (DLL4), ELISA Kit

Cat PEK02373

Delta-like protein 4

UniProt Synonym Protein Names

Drosophila Delta homolog 4; Delta4

UniProt Entry Name

DLL4_HUMAN

UniProt Secondary Accession

Q3KP23; Q9NQT9

UniProt Related Accession

Q9NR61

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY