

Sucrose Phosphate Synthase (SPS), ELISA Kit

Cat PEK02387

Species

Plant

Full Product Name

Plant Sucrose Phosphate Synthase (SPS) ELISA Kit

Product Gene Name

SPS elisa kit

Assay Type

Sandwich

Detection Range

15.6pmol/L-500pmol/L.

Sensitivity

2.0pmol/L.

Intra-assay Precision

Intra-assay CV (%) is less than 15%.

Inter-assay Precision

Inter-assay CV (%) is less than 15%. [CV(%) = SD/mean ×100].

Specificity

No significant cross-reactivity or interference between this analyte and analogues is observed.

Samples

Plant Tissue Homogenate and other Plant Biological Fluids samples

Storage

Store all reagents at 2-8°C

Protein Family

Probable sucrose-phosphate synthase

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY