

# BH3 Interacting Domain Death Agonist (BID), ELISA Kit

Cat PEK02598

---

## Species

Plant

## Full Product Name

Plant BH3 Interacting Domain Death Agonist (BID) ELISA Kit

## Product Gene Name

BID elisa kit

## Molecular Weight

11,263 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

BH3-interacting domain death agonist

## NCBI Accession #

AAC34365.1

## NCBI GI #

3540247

## NCBI GeneID

637

## NCBI Official Full Name

BH3 interacting domain death agonist

## NCBI Official Symbol

BID

## NCBI Official Synonym Full Names

BH3 interacting domain death agonist

## NCBI Official Synonym Symbols

FP497

## NCBI Protein Information

BH3-interacting domain death agonist

## UniProt Gene Name

BID

## UniProt Synonym Gene Names

BID

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

# BH3 Interacting Domain Death Agonist (BID), ELISA Kit

Cat PEK02598

---

BH3-interacting domain death agonist

## UniProt Synonym Protein Names

p22 BID; BID

## UniProt Entry Name

BID\_HUMAN

## UniProt Secondary Accession #

Q549M7; Q71T04; Q7Z4M9; Q8IY86

## UniProt Related Accession #

P55957

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**