

Signal Transducer and Activator of Transcription 5B (STAT5B), ELISA Kit

Cat PEK02607

Species

Plant

Full Product Name

Plant Signal Transducer and Activator of Transcription 5B (STAT5B) ELISA Kit

Product Gene Name

STAT5B elisa kit

Molecular Weight

89,866 Da

Storage

Store all reagents at 2-8°C

Protein Family

Signal transducer and activator of transcription

NCBI Accession

NP_036580.2

NCBI GI

21618344

NCBI GeneID

6777

NCBI Official Full Name

signal transducer and activator of transcription 5B

NCBI Official Symbol

STAT5B

NCBI Official Synonym Full Names

signal transducer and activator of transcription 5B

NCBI Official Synonym Symbols

STAT5

NCBI Protein Information

signal transducer and activator of transcription 5B

UniProt Gene Name

STAT5B

UniProt Protein Name

Signal transducer and activator of transcription 5B

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Signal Transducer and Activator of Transcription 5B (STAT5B), ELISA Kit

Cat PEK02607

UniProt Entry Name

STA5B_HUMAN

UniProt Secondary Accession

Q8WWS8

UniProt Related Accession

P51692

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY