

Major Vault Protein (MVP), ELISA Kit

Cat PEK02618

Species

Plant

Full Product Name

Plant Major Vault Protein (MVP) ELISA Kit

Product Gene Name

MVP elisa kit

Molecular Weight

95,670 Da

Storage

Store all reagents at 2-8°C

Protein Family

Hyperpolarization-activated voltage-gated potassium channel

NCBI Accession

NP_001079939.1

NCBI GI

147903423

NCBI GeneID

379630

NCBI Official Full Name

major vault protein

NCBI Official Symbol

mvp

NCBI Official Synonym Full Names

major vault protein

NCBI Protein Information

major vault protein

UniProt Gene Name

mvp

UniProt Synonym Gene Names

MVP

UniProt Protein Name

Major vault protein

UniProt Entry Name

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY

Major Vault Protein (MVP), ELISA Kit

Cat PEK02618

MVP_XENLA

UniProt Related Accession

Q6PF69

FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY