

# Leukotriene B4 Receptor 2 (LTB4R2), ELISA Kit

Cat PEK02622

---

## Species

Plant

## Full Product Name

Plant Leukotriene B4 Receptor 2 (LTB4R2) ELISA Kit

## Product Gene Name

LTB4R2 elisa kit

## Molecular Weight

37,942 Da

## Storage

Store all reagents at 2-8°C

## Protein Family

Leukotriene B4 receptor

## NCBI Accession #

NP\_062813.2

## NCBI GI #

291042665

## NCBI GeneID

56413

## NCBI Official Full Name

leukotriene B4 receptor 2

## NCBI Official Symbol

LTB4R2

## NCBI Official Synonym Full Names

leukotriene B4 receptor 2

## NCBI Official Synonym Symbols

BLT2; NOP9; BLTR2; JULF2; KPG\_004; LTB4-R2; LTB4-R 2

## NCBI Protein Information

leukotriene B4 receptor 2

## UniProt Gene Name

LTB4R2

## UniProt Synonym Gene Names

BLT2R; BLTR2; LTB4-R 2; LTB4-R2

## UniProt Protein Name

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**

---

# Leukotriene B4 Receptor 2 (LTB4R2), ELISA Kit

Cat PEK02622

---

Leukotriene B4 receptor 2

## UniProt Synonym Protein Names

LTB4 receptor JULF2; Leukotriene B4 receptor BLT2; Seven transmembrane receptor BLTR2

## UniProt Entry Name

LT4R2\_HUMAN

## UniProt Secondary Accession #

Q5KU28; Q9NPE5

## UniProt Related Accession #

Q9NPC1

---

**FOR RESEARCH OR FURTHER MANUFACTURING USE ONLY**